

	Performance Hyperbarics	Competition Chamber large	Competition Chamber medium	Competition Chamber small
Length	105 inches	92	90	76 inches
Diameter	34 inches	32	27	22 inches
Volume	95,283 inches ³	73,953 inches ³	51,504 inches ³	28,875 inches ³
Zipper length	65 inches	55 inches	55 inches	55 inches
Windows	4	4	3	3
Additional Accessory Ports	4	1	1	1
Internal Metal Frame Outer Diameter	1 1/8th inch	1/4 inch	1/4 inch	1/4 inch
Factory test level	8 p.s.i.	6 p.s.i.	6 p.s.i.	6 p.s.i.
Solo operation possible	Yes	Yes	Yes	Yes
Zipper material	Metal	Plastic	Plastic	Plastic
Anti-roll bolsters	Yes	Yes	Yes	Yes
Oxygen Safe Intercom Option	Yes	No	No	No
Double Bag design	Yes	No	No	No
Reinforcement Straps and Rings	Yes	No	No	No
Oxygen Analysis Port	Yes	No	No	No
Internal Pressure Gauge Option	Yes	No	No	No
Compressor cooler and dehumidifier	Yes	No	No	No
Quiet air input system	Yes	No	No	No
Price	\$14,900	\$20,900	\$16,900	\$10,500